#### **Blueberry Project**

#### Sandra Nnadi

The Bluebe Plant Introduction

Characterization

# Understanding The Fungal Microbiome in Blueberry Roots

#### Sandra Nnadi

University of Vermont Department of Plant Biology Burlington, VT 05401 USA

sandra.nnadi@uvm.edu

3 March 2021



### **Overview**

**Blueberry Project** 

Sandra Nnadi

The Blueber Plant Introduction

Molecular Characterizatio

- 1 The Blueberry Plant
  - Introduction

2 Molecular Characterization

## The Blueberry Plant

**Blueberry Project** 

Sandra Nnadi

The Blueberry Plant

Introduction

Molecular Characterization

- Northern Highbush Blueberry
- Ericaceae family
- grows in acidic soil
- mycorrhizae may improve performance
- composition of fungal community is poorly understood

# **Figure**

#### **Blueberry Project**

Sandra Nnadi

The Blueber Plant Introduction

Molecular Characterization



### **Table**

**Blueberry Project** 

Sandra Nnadi

The Blueberry

Introduction

Molecular Characterizati

Treatments	Inoculation	Fertilization
Control	No	No
Commercial Inoculum	Yes	No
Native Inoculum	Yes	Yes

Table: Experimental Design

### **Molecular Characterization**

Blueberry Project

Sandra Nnadi

The Blueberry Plant

Molecular Characterization A Molecular approach will be used to understand the composition of the fungal community. The Internal Transcribed Spacer ITS is located between the small and large subunit of the ribosomal RNA.

The ribosomal RNA is highly conserved in all fungal species but the ITS region varies and this variation will be used to distinguish species within samples.

# Library Preparation

1st PCR and Cleanup

Blueberry Project

Sandra Nnadi

The Blueberry

Introduction

Molecular

Characterization

beads.

2nd PCR and cleanup

Add Index primers to adapter overhang to provide a unique identifier for each sample before pooling then clean with Ampure XP beads.

Amplify the ITS1 region with ITS forward and reverse

primers then clean the PCR product with Ampure XP

**Pooling and Sequencing** 

Check quality of samples with Bioanalyzer trace, dilute final Library to 4nM then pool samples in a 2ml screw cap and submit for Illumina MiSeg Seguencing.

#### **Blueberry Project**

#### Sandra Nnadi

The Blueberry
Plant
Introduction

Molecular Characterization

# The End